



WHERE YOU WANT TO BE.

Sapphire[®] 2





Sapphire[®] 2

Take it apart,
take it away,
take it further!





Detachable

To get your Sapphire 2 ready for transporting, simply lift off the seat, lift out the batteries, and lift the rear handle to separate the front and rear of the scooter. It really is that simple! (and by the way it goes back together just as smoothly!).



The Sterling Lock system is cable and connector free so there are no fiddly bits to deal with.

Transportable

Now the Sapphire 2 can be loaded into your vehicle. Luckily, as one of the lightest scooters in its class, this is no problem with the Sterling Lift dual handles fitted to the batteries and rear drive unit. These make lifting and handling the parts safer and more comfortable.



With its large batteries carrying you in style and comfort, the Sapphire 2 is the scooter you'll want to go further on!


The tiller can be safely secured in position for transporting.

Comfortable

Now you're ready! Stretch out your legs and take advantage of the Sterling Comfort package: three wheel legroom with four wheel stability, infinitely adjustable tiller, wrap-around handle bars and finger-thumb directional control all fitted as standard.



Sapphire[®] 2



Technical data sheet

Version	4 Wheels	Drive unit weight	19.5kg
Speed	4mph (6.4kph)	Battery weight	unit 12kg
Batteries	35 Amph	Total weight when folded/dismantled	22.5kg
Charger off board	5 Amph	Turning radius	140cm (56")
Maximum range	30km (19 miles)	Maximum safe slope	14° - 113 kg user weight / 10° - 150kg user weight
Colour	Blue & Red	Ground clearance	9.5cm (3.5")
Max user weight	150kg (23 stone)	Tyres	Pneumatic
Total weight	79kg	Front wheel size	22.5cm (9")
Heaviest part (front section)	22.5kg	Rear wheel size	26.5cm (10.5")
Seat weight	12.5kg	Stylish lights	Yes

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



PB_SAPPHIRE2_EN_GB_REV1

Sunrise Medical Limited
High Street, Wollaston
Stourbridge
West Midlands
DY8 4PS

Tel: +44 (0) 1384 44 66 66
Fax: +44 (0) 1384 44 66 44

www.SunriseMedical.com