



### Get more than mobility

The *Invacare TDX SP* is a high-performance power chair for those that want to get the maximum out of their wheelchair. The *TDX SP* features Centre Wheel Drive technology for outstanding manoeuvrability. Active customers will appreciate the Stability Lock function which offers maximum stability and safety on uneven terrain and SureStep®, for smooth travel over obstacles. Our pioneering tilt-in-space function provides optimum pressure relief and positioning. Even heavier customers can enjoy speedy, secure mobility thanks to a powerful drive mode.



**Yes, you can.®**



**Yes to mobility, no to compromise**

The tightest possible turning radius is based in the Centre Wheel Drive construction. This type of wheel base provides excellent manoeuvrability and an intuitive driving experience. The *Invacare TDX SP* is the ideal companion for both indoor and outdoor conditions thanks to good stability and high traction on uneven ground.

**Tried, tested .... and tough**

For maximum stability and safety, the patented Stability Lock function keeps all six wheels on the ground. The SureStep® feature enables a smooth ride over even the toughest surfaces. For transportation, the crash-tested 'tie-down' ensures the chair is transported securely.

**Designed for heavier customers**

Thanks to sophisticated engineering, robust construction and powerful motors, the *TDX SP* is ideal for heavy-duty use and can support a weight of up to 180 kg.

**Features and Options**



**Centre Wheel Drive**

This specific driving technology allows excellent manoeuvrability in narrow spaces and offers a responsive driving experience.



**SureStep**

Original patented technology carries the wheelchair up and over obstacles up to 7.5 cm high, and guarantees a smooth downwards transition.



**Stability Lock**

For maximum stability and safety, this innovative feature keeps all six wheels on the ground and also has a shock-absorbing effect.

**Features and Options**



**Seat lifter**  
The **TDX SP** features a seat lifter of 30 cm including powered tilt function of up to 28°.



**Tilt-in-space**  
A maximum seat angle of 28° ensures a comfortable and safe sitting position.



**Reclining back rest**  
A maximum backrest angle of 35° provides excellent seating and back support.



**4-point tie-down**  
Four strong hooks ensure a stable position during transportation. The **TDX SP** has been successfully crash-tested.

## Technical data

Standard Seat	44 - 48 cm	41 / 46 / 51 cm	45 cm <sup>1</sup> / 49 cm <sup>7</sup>	48 / 54 cm <sup>2</sup>	23 - 32 cm	29 - 47 cm <sup>3</sup>	0 - 28° <sup>7</sup>			
Contour Seat	44 - 53 cm	41 / 46 / 51 cm	45 cm <sup>1</sup> / 49 cm <sup>7</sup>	64 cm <sup>2</sup>	23 - 32 cm	29 - 47 cm <sup>3</sup>	0 - 28° <sup>7</sup>			
Flex Seat	44 - 53 cm	38 / 43 / 48 cm	45 cm <sup>1</sup> / 49 cm <sup>7</sup>	55 cm <sup>2</sup>	23 - 32 cm	29 - 47 cm <sup>3</sup>	0 - 28° <sup>7</sup>			
Standard Seat	-10° - 30°	64 cm <sup>4</sup>	98 cm <sup>5</sup>	98 cm <sup>6</sup>	91 cm	138 kg <sup>5</sup>	180 kg / 160 kg <sup>7</sup> / 140 kg <sup>8</sup>			
Contour Seat	0° - 45°	64 cm <sup>4</sup>	108 cm <sup>5</sup>	98 cm <sup>6</sup>	91 cm	138 kg <sup>5</sup>	180 kg / 160 kg <sup>7</sup> / 140 kg <sup>8</sup>			
Flex Seat	0° - 45° <sup>9</sup>	64 cm <sup>4</sup>	108 cm <sup>5</sup>	98 cm <sup>6</sup>	91 cm	138 kg <sup>5</sup>	180 kg / 160 kg <sup>7</sup> / 140 kg <sup>8</sup>			
	73,6 Ah (C20)	112 cm	7,5 cm	9° / 16 %	2 x 340 W	39 km <sup>10</sup>	8 km/h and 6 km/h			



Bariatric use



G-Trac™ option available

- 1/ ground to seat plate without cushion
- 2/ measured from seat plate
- 3/ measured from foot plate to seat plate
- 4/ depending on seat adjustment
- 5/ standard configuration, varies with options

- 6/ leg rests at largest angle
- 7/ with seat lifter or powered tilt
- 8/ programmed to 6 km/h
- 9/ at seat depth of 48 cm the powered recliner is limited to 30°
- 10/ according to ISO 7176-4:2008

## Frame colour



Electric Red

INVACARE® Limited  
 Pencoed Technology Park - Pencoed -  
 Bridgend CF35 5HZ  
 Tel: +44 1 656 776 222  
 Fax: +44 1656 776 220  
 E-mail: uk@invacare.com  
 www.invacare.co.uk  
 Sales Order E-mail:  
 ordersuk@invacare.com

INVACARE® Ireland Ltd  
 Unit 5 - Seatown Business Campus  
 Seatown Road - SWORDS - County  
 Dublin  
 Tel. +353 1 8107084  
 Fax +353 1 8107085  
 Email: Ireland@invacare.com  
 www.invacare.ie

© 2010 Invacare® International Sarl  
 All rights reserved.

All information quoted is believed to be  
 correct at time of print.  
 Invacare® reserves the right to alter  
 product specifications without prior  
 consultation.

TDX SP - UK - 11/2010 - 150762



More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.



Yes, you can.®